

Product datasheet for **HP233467**

FCGR1CP Human qPCR Primer Pair (NM_001128589)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	100132417
Forward Sequence:	ATACAGGTGCCAGAGAGGTCTC
Reverse Sequence:	CCAGCTTATCCTTCCACGCATG
ACCN:	NM_001128589 , NM_001128589.1 , BF061349
Synonyms:	CD64; Fc fragment of IgG, high affinity Ic, receptor (CD64); Fc fragment of IgG, high affinity Ic, receptor for (CD64); IGFRC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)