

Product datasheet for **HP233499**

IL15 Human qPCR Primer Pair (NM_172174)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3600
Forward Sequence:	AACAGAAGCCAACCTGGGTGAATG
Reverse Sequence:	CTCCAAGAGAAAGCACTTCATTGC
ACCN:	<u>NM_172174</u> , <u>NM_172174.1</u> , <u>NM_172174.2</u> , <u>BC018149</u> , <u>BC018149.2</u> , <u>BC100961</u> , <u>BC100962</u> , <u>BC100963</u> , <u>BX116580</u>
UniProt ID:	<u>P40933</u>
Synonyms:	IL-15; Interleukin 15; MGC9721; OTTHUMP00000164617
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)