

Product datasheet for **HP233547**

Caveolin 3 (CAV3) Human qPCR Primer Pair (NM_033337)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 859 |
| Forward Sequence: | GATTGACCTGGTGAACCGAGAC |
| Reverse Sequence: | CTTGAGACAGTGAAGGTGGTG |
| ACCN: | NM_033337 , NM_033337.1 , NM_033337.2 , BC069368 , BC102036 , BC102033 , BC102037 , BF310794 , BX114193 , NM_033337.3 |
| UniProt ID: | P56539 |
| Synonyms: | LGMD1C; LQT9; MPDT; RMD2; VIP-21; VIP21 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)