

## Product datasheet for **HP233617**

### CCDC169 Human qPCR Primer Pair (NM\_001144981)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	728591
Forward Sequence:	GAACGCCGTACATACCTAGCTG
Reverse Sequence:	CTCACTGTCTTCTGCTTAGCTGG
ACCN:	<u><a href="#">NM_001144981</a></u> , <u><a href="#">NM_001144981.1</a></u> , <u><a href="#">NM_001144981.2</a></u> , <u><a href="#">BC146975</a></u> , <u><a href="#">BC037403</a></u> , <u><a href="#">BC171867</a></u> , <u><a href="#">BC171872</a></u> , <u><a href="#">BC171874</a></u> , <u><a href="#">BC171880</a></u> , <u><a href="#">BC171881</a></u> , <u><a href="#">BC171882</a></u> , <u><a href="#">BC171884</a></u> , <u><a href="#">NM_001144981.3</a></u>
UniProt ID:	<u><a href="#">A6NNP5</a></u>
Synonyms:	C13orf38
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)