

Product datasheet for **HP233622**

ATP6V0A4 Human qPCR Primer Pair (NM_130840)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	50617
Forward Sequence:	TTGGAGAGCGTCAATGTGAGGC
Reverse Sequence:	CTTCTGCACCTTGATGAGCCAG
ACCN:	<u>NM_130840</u> , <u>NM_130840.1</u> , <u>NM_130840.2</u> , <u>BC109304</u> , <u>BC109305</u> , <u>NM_130840.3</u>
UniProt ID:	<u>Q9HBG4</u>
Synonyms:	A4; ATP6N1B; ATP6N2; DRTA3; RDRTA2; RTA1C; RTADR; STV1; VPH1; VPP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)