

Product datasheet for **HP233654**

NUP43 Human qPCR Primer Pair (NM_198887)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	348995
Forward Sequence:	GACGTTGCTACAGGATCTTG
Reverse Sequence:	CCTCAAACCTCCATCAGAGTC
ACCN:	<u>NM_198887</u> , <u>NM_198887.1</u> , <u>NM_198887.2</u> , <u>BC065028</u> , <u>BC065028.1</u> , <u>BC040121</u> , <u>BC046637</u> , <u>BC047539</u> , <u>BG527619</u> , <u>BG567254</u> , <u>BX647097</u> , <u>NM_198887.3</u>
UniProt ID:	<u>Q8NFH3</u>
Synonyms:	bA350J20.1; p42
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)