

## Product datasheet for **HP233665**

### **IQCF5 Human qPCR Primer Pair (NM\_001145059)**

#### **Product data:**

|                          |                                                                                                                                                                                                                                                                                                                                              |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b>     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                            |
| <b>Gene ID:</b>          | 389124                                                                                                                                                                                                                                                                                                                                       |
| <b>Forward Sequence:</b> | TGGAGAAGCTGCTGGCAAAGAG                                                                                                                                                                                                                                                                                                                       |
| <b>Reverse Sequence:</b> | GAGCAAACGACAGTAACGCTGG                                                                                                                                                                                                                                                                                                                       |
| <b>ACCN:</b>             | <u><a href="#">NM_001145059</a></u> , <u><a href="#">NM_001145059.1</a></u> , <u><a href="#">BX106781</a></u> , <u><a href="#">NM_001145059.2</a></u>                                                                                                                                                                                        |
| <b>UniProt ID:</b>       | <u><a href="#">A8MTL0</a></u>                                                                                                                                                                                                                                                                                                                |
| <b>Components:</b>       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)                                                                                                                                                                                                                                               |
| <b>Quality Control:</b>  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| <b>Storage:</b>          | Store at -20°C.                                                                                                                                                                                                                                                                                                                              |
| <b>Stability:</b>        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                 |



[View online »](#)