

## Product datasheet for **HP233722**

### **AZ11 (CEP131) Human qPCR Primer Pair (NM\_014984)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	22994
<b>Forward Sequence:</b>	GTCCTGAGTGAAAAGTGCAGG
<b>Reverse Sequence:</b>	GCTTTCTCGGTGGCGCTCATTA
<b>ACCN:</b>	<a href="#">NM_014984</a> , <a href="#">NM_014984.1</a> , <a href="#">NM_014984.2</a> , <a href="#">NM_014984.3</a> , <a href="#">BC010413</a> , <a href="#">BC011615</a> , <a href="#">BC012130</a> , <a href="#">BC040023</a> , <a href="#">BC136578</a> , <a href="#">BC136580</a> , <a href="#">BC144263</a> , <a href="#">BE732160</a> , <a href="#">BE734213</a> , <a href="#">BG773915</a> , <a href="#">BU742835</a> , <a href="#">BX394713</a>
<b>UniProt ID:</b>	<a href="#">Q9UPN4</a>
<b>Synonyms:</b>	AZ1; AZ11; ZA1
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)