

Product datasheet for **HP233723**

C1orf124 (SPRTN) Human qPCR Primer Pair (NM_032018)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83932
Forward Sequence:	GTGTTCCATCCGTCTCAGCGAA
Reverse Sequence:	ATGCCCTTCTCGGTCTTTGTCG
ACCN:	NM_032018 , NM_032018.1 , NM_032018.2 , NM_032018.3 , NM_032018.4 , NM_032018.5 , NM_032018.6 , BC068478 , BC068478.1 , BC015740 , BM480135 , BU630543
UniProt ID:	Q9H040
Synonyms:	C1orf124; DVC1; PRO4323; spartan
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)