

Product datasheet for **HP233730**

MRP4 (ABCC4) Human qPCR Primer Pair (NM_005845)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10257
Forward Sequence:	CTGTTGGAGGATGGTGATCTGAC
Reverse Sequence:	CTGCTAACTTCCGCATCTACTGC
ACCN:	NM_005845 , NM_005845.1 , NM_005845.2 , NM_005845.3 , NM_005845.4 , BC041560 , BM687848 , BP382967 , BU689822 , BX113883 , BX647118
UniProt ID:	O15439
Synonyms:	MOAT-B; MOATB; MRP4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)