

Product datasheet for **HP233735**

Piccolo (PCLO) Human qPCR Primer Pair (NM_033026)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27445
Forward Sequence:	TCTGACTTCAGGGAGAGTTACTAC
Reverse Sequence:	GGTGGTGCATTCTTGCTGTGG
ACCN:	NM_033026 , NM_033026.1 , NM_033026.2 , NM_033026.3 , NM_033026.4 , NM_033026.5 , BC001304 , BC108289 , BC122565 , BC125271 , BE538315 , BG913494 , BX649041 , NM_033026.6
UniProt ID:	Q9Y6V0
Synonyms:	ACZ; PCH3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)