

Product datasheet for **HP233741**

RNF34 Human qPCR Primer Pair (NM_194271)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80196
Forward Sequence:	GGAAGCTGGTAGAGAAAGTAAACCG
Reverse Sequence:	ACACAGGCTGTCGTCTTCTCA
ACCN:	<u>NM_194271</u> , <u>NM_194271.1</u> , <u>NM_194271.2</u> , <u>BC007826</u> , <u>BC007826.2</u> , <u>BC039029</u> , <u>BC046174</u> , <u>BI559962</u> , <u>BT007283</u>
UniProt ID:	<u>Q969K3</u>
Synonyms:	CARP-1; CARP1; hRFI; RFI; RIF; RIFF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)