

Product datasheet for **HP233776**

B3GALT5 Human qPCR Primer Pair (NM_033172)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10317
Forward Sequence:	AGCGGAAACGAAAGAGGTGGAC
Reverse Sequence:	CCTGAGGACAAAAGCGATGGAC
ACCN:	<u>NM_033172</u> , <u>NM_033172.1</u> , <u>NM_033172.2</u> , <u>BC104864</u> , <u>BC104864.1</u> , <u>BC104862</u>
UniProt ID:	<u>Q9Y2C3</u>
Synonyms:	3-GalTase 5; B3GalT-V; B3GalTx; B3T5; beta-1; beta-3-Gx-T5; beta3Gal-T5; GLCT5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)