

Product datasheet for **HP233820**

PCDHA6 Human qPCR Primer Pair (NM_031848)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56142
Forward Sequence:	CTGCCACATCTTCACGGTGTCT
Reverse Sequence:	GCACCGAAATGTAGCTCGACAAC
ACCN:	<u>NM_031848</u> , <u>NM_031848.1</u> , <u>NM_031848.2</u> , <u>BC036674</u> , <u>BC036674.2</u> , <u>BC019862</u> , <u>BC042887</u> , <u>BC043441</u> , <u>NM_031848.3</u>
UniProt ID:	<u>Q9UN73</u>
Synonyms:	CNR2; CNRN2; CNRS2; CRNR2; PCDH-ALPHA6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)