

Product datasheet for **HP233857**

TNRC6C Human qPCR Primer Pair (NM_001142640)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57690
Forward Sequence:	CACCGTTGCTTGGTCCAGTTTC
Reverse Sequence:	TGCTGTGCCATTATCCACCAGG
ACCN:	NM_001142640 , NM_001142640.1 , BC035105 , BC039479 , BC045631 , BC051468 , BC156366 , BC157115
UniProt ID:	Q9HCJ0
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)