

Product datasheet for **HP233858**

MS4A3 Human qPCR Primer Pair (NM_006138)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	932
Forward Sequence:	GGGTGTCTTTCTGGGTCCTTG
Reverse Sequence:	CCTGCTACAACAGACAAGTTCC
ACCN:	NM_006138 , NM_006138.1 , NM_006138.2 , NM_006138.3 , NM_006138.4 , BC008487 , BC008487.1 , BE564041 , BF210971 , BG495577 , NM_006138.5
UniProt ID:	Q96HJ5
Synonyms:	CD20L; HTM4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)