

Product datasheet for **HP233883**

Siglec 7 (SIGLEC7) Human qPCR Primer Pair (NM_016543)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27036
Forward Sequence:	CTGGTCTTCCTCTCCTTCTGTG
Reverse Sequence:	GCATCCTTCATGCCTATGTCTCC
ACCN:	<u>NM_016543</u> , <u>NM_016543.1</u> , <u>NM_016543.2</u> , <u>NM_016543.3</u> , <u>BC028150</u> , <u>BC028150.1</u> , <u>BX644517</u>
UniProt ID:	<u>Q9Y286</u>
Synonyms:	AIRM-1; AIRM1; CD328; CDw328; D-siglec; p75; p75/AIRM1; QA79; SIGLEC-7; SIGLEC19P; SIGLECP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

