

## Product datasheet for **HP233888**

### A1CF Human qPCR Primer Pair (NM\_138932)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29974
Forward Sequence:	CAGCCACCAAAGGACATCTCAG
Reverse Sequence:	CTGTGTATGCCAAATAGCCACGG
ACCN:	<a href="#">NM_138932</a> , <a href="#">NM_138932.1</a> , <a href="#">NM_138932.2</a> , <a href="#">BC130519</a> , <a href="#">BC022263</a> , <a href="#">BC054873</a> , <a href="#">BC144194</a> , <a href="#">BC144196</a> , <a href="#">BG188383</a>
UniProt ID:	<a href="#">Q9NQ94</a>
Synonyms:	ACF; ACF64; ACF65; APOBEC1CF; ASP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)