

Product datasheet for **HP233897**

GTPBP10 Human qPCR Primer Pair (NM_033107)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	85865
Forward Sequence:	GTGGTGATGTCTGGGTTGTAGC
Reverse Sequence:	CACAGTCTTTTCCTTTGGAGCCT
ACCN:	NM_033107 , NM_033107.1 , NM_033107.2 , NM_033107.3 , BC107714 , BC107714.1 , BC004923 , BC021573 , BC144162 , NM_033107.4
UniProt ID:	A4D1E9
Synonyms:	ObgH2; UG0751c10
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)