

Product datasheet for **HP233921**

Pregnancy specific beta 1 glycoprotein 11 (PSG11) Human qPCR Primer Pair (NM_002785)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5680
Forward Sequence:	AAGCGAGGTGATGGGACTAGAG
Reverse Sequence:	GGTTAAGATCACAGTCTCCATGG
ACCN:	NM_002785 , NM_002785.1 , NM_002785.2 , BC020711 , BC020711.1 , BE767872 , BF084984 , BG571293 , BG619770 , NM_002785.3
UniProt ID:	Q9UQ72
Synonyms:	PSBG-11; PSBG-13; PSG13; PSG14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)