

## Product datasheet for **HP233954**

### FMNL2 Human qPCR Primer Pair (NM\_052905)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	114793
Forward Sequence:	TGCCTGGAACATCTTCACCCAC
Reverse Sequence:	CTGATTGGGCTTCAGAGCAACC
ACCN:	<u><a href="#">NM_052905</a></u> , <u><a href="#">NM_052905.1</a></u> , <u><a href="#">NM_052905.2</a></u> , <u><a href="#">NM_052905.3</a></u> , <u><a href="#">BC014343</a></u> , <u><a href="#">BC036492</a></u> , <u><a href="#">BC113878</a></u> , <u><a href="#">BC114438</a></u> , <u><a href="#">BC167159</a></u> , <u><a href="#">BC167804</a></u> , <u><a href="#">BQ067763</a></u> , <u><a href="#">BU633028</a></u> , <u><a href="#">NM_052905.4</a></u>
UniProt ID:	<u><a href="#">Q96PY5</a></u>
Synonyms:	FHOD2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)