

## Product datasheet for **HP233960**

### U2AF35 (U2AF1) Human qPCR Primer Pair (NM\_006758)

#### Product data:

|                   |   |
|-------------------|---|
| Product Type:     | qPCR Primer Pairs   |
| Gene ID:          | 7307  |
| Forward Sequence: | GAACGTCTGTGACAACCTGGGA  |
| Reverse Sequence: | AGGCTTCTCTGAAGTCCGTAC   |
| ACCN:             | <u><a href="#">NM_006758</a></u> , <u><a href="#">NM_006758.1</a></u> , <u><a href="#">NM_006758.2</a></u> , <u><a href="#">BC001177</a></u> , <u><a href="#">BC005915</a></u> , <u><a href="#">BC005915.1</a></u> , <u><a href="#">BC001923</a></u> , <u><a href="#">BC013786</a></u> , <u><a href="#">BC032521</a></u> , <u><a href="#">NM_006758.3</a></u> |
| UniProt ID:       | <u><a href="#">Q01081</a></u>   |
| Synonyms:         | FP793; RN; RNU2AF1; U2AF35; U2AFBP  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.                  |
| Storage:          | Store at -20°C.   |
| Stability:        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)