

Product datasheet for **HP233964**

Integrin alpha 3 (ITGA3) Human qPCR Primer Pair (NM_005501)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3675
Forward Sequence:	GCCTGACAACAAGTGTGAGAGC
Reverse Sequence:	GGTGTTTCGTCACGTTGATGCTC
ACCN:	<u>NM_005501</u> , <u>NM_005501.1</u> , <u>NM_005501.2</u> , <u>BC015344</u> , <u>BC043917</u> , <u>BC136636</u> , <u>BC144328</u> , <u>BC150190</u> , <u>BM978115</u> , <u>BU631476</u>
UniProt ID:	<u>P26006</u>
Synonyms:	CD49C; GAP-B3; GAPB3; ILNEB; MSK18; VCA-2; VL3A; VLA3a
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)