

## Product datasheet for **HP233971**

### RBMS3 Human qPCR Primer Pair (NM\_001003793)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27303
Forward Sequence:	GTGTTGGCTTTGCCAGAATGGAG
Reverse Sequence:	CGCTTCTTTTGCCTCCATCAGC
ACCN:	<u><a href="#">NM_001003793</a>, <a href="#">NM_001003793.1</a>, <a href="#">NM_001003793.2</a>, <a href="#">BC117315</a>, <a href="#">BC143521</a>, <a href="#">NM_001003793.3</a></u>
UniProt ID:	<u><a href="#">Q6XE24</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)