

Product datasheet for **HP233992**

SRRM3 Human qPCR Primer Pair (NM_001110199)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	222183
Forward Sequence:	AGACTGGAGTCCGAATGCAGCT
Reverse Sequence:	ATCTGTGCCGTCTTCTCCTCCG
ACCN:	<u>NM_001110199</u> , <u>NM_001110199.1</u> , <u>NM_001110199.2</u> , <u>BC063788</u> , <u>BC029292</u> , <u>BC045627</u> , <u>BC071856</u> , <u>BC096755</u> , <u>BE466083</u> , <u>NM_001110199.3</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)