

## Product datasheet for **HP234007**

### HMBOX1 Human qPCR Primer Pair (NM\_001135726)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79618
Forward Sequence:	GCAGCAATCCTGGAGAGTCATG
Reverse Sequence:	GTCTTGGTGGTCACTATGGCTC
ACCN:	<a href="#">NM_001135726</a> , <a href="#">NM_001135726.1</a> , <a href="#">NM_001135726.2</a> , <a href="#">BC069242</a> , <a href="#">BC069242.1</a> , <a href="#">BC009259</a> , <a href="#">BF369137</a> , <a href="#">BQ019252</a> , <a href="#">BX353564</a> , <a href="#">BX648950</a> , <a href="#">NM_001135726.3</a>
UniProt ID:	<a href="#">Q6NT76</a>
Synonyms:	HNF1LA; HOT1; PBHNF; TAH1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)