

Product datasheet for **HP234056**

CD239 (BCAM) Human qPCR Primer Pair (NM_005581)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4059
Forward Sequence:	GATGGCATCTCCTGTGAAGCCT
Reverse Sequence:	CACGCAGTAGAAGACAGCAACG
ACCN:	<u>NM_005581</u> , <u>NM_005581.2</u> , <u>NM_005581.3</u> , <u>NM_005581.4</u> , <u>BC050450</u> , <u>BC050450.1</u> , <u>BC036745</u> , <u>BC128254</u> , <u>BM985314</u> , <u>BP352029</u> , <u>BU838291</u> , <u>BX387680</u> , <u>NM_005581.5</u>
UniProt ID:	<u>P50895</u>
Synonyms:	AU; CD239; LU; MSK19
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

