

Product datasheet for **HP234110**

SMIM29 Human qPCR Primer Pair (NM_178508)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	221491
Forward Sequence:	AGCTGAGGAAGTGCATGAGGCT
Reverse Sequence:	AGGTCACGTCTTGACATCCAGC
ACCN:	<u>NM_178508</u> , <u>NM_178508.1</u> , <u>NM_178508.2</u> , <u>NM_178508.3</u> , <u>NM_178508.4</u> , <u>BC023627</u> , <u>BC023627.2</u> , <u>BC047919</u> , <u>BC047919.1</u> , <u>BC009526</u> , <u>BM474526</u> , <u>BU194178</u> , <u>BU956518</u>
UniProt ID:	<u>Q86T20</u>
Synonyms:	C6orf1; LBH
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)