

Product datasheet for **HP234147**

CHST15 Human qPCR Primer Pair (NM_015892)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51363
Forward Sequence:	ATCAAGGAGCCCACTGGTGGGA
Reverse Sequence:	GCAGTCCTTGATGGATCTGGTG
ACCN:	NM_015892 , NM_015892.1 , NM_015892.2 , NM_015892.3 , NM_015892.4 , BC027908 , BC027908.1 , BC050540 , BC075813 , BM193389 , BM977716 , NM_015892.5
UniProt ID:	Q7LFX5
Synonyms:	BRAG; GALNAC4S-6ST
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)