

Product datasheet for **HP234176**

GFI1 Human qPCR Primer Pair (NM_001127216)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2672
Forward Sequence:	GCTTCAAGAGGTCATCCACACTG
Reverse Sequence:	ACCTGGCACTTGTGAGGCTTCT
ACCN:	NM_001127216 , NM_001127216.1 , BC032751 , BC032751.1 , BC074866 , BC074867 , BE242676 , BM150826 , BQ336135 , NM_001127216.2
UniProt ID:	Q99684
Synonyms:	GFI-1; GFI1A; SCN2; ZNF163
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)