

Product datasheet for **HP234194**

FTSJD1 (CMTR2) Human qPCR Primer Pair (NM_018348)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55783
Forward Sequence:	CTGACCACTCTTGAAACGGTG
Reverse Sequence:	GACTTCGGAGTTTCCTGCCTTG
ACCN:	<u>NM_018348</u> , <u>NM_018348.1</u> , <u>NM_018348.2</u> , <u>NM_018348.3</u> , <u>NM_018348.4</u> , <u>NM_018348.5</u> , <u>BC035005</u> , <u>BC035005.2</u> , <u>BC034468</u> , <u>BU429551</u> , <u>BX648660</u> , <u>NM_018348.6</u>
UniProt ID:	<u>Q8IYT2</u>
Synonyms:	AFT; FTSJD1; HMTr2; MTr2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)