

Product datasheet for **HP234228**

Caldesmon (CALD1) Human qPCR Primer Pair (NM_033138)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 800 |
| Forward Sequence: | CTGTTCTGCTGAAGGTGTACG |
| Reverse Sequence: | CCTACCTCAAGCCAGCAGTTTC |
| ACCN: | NM_033138 , NM_033138.1 , NM_033138.2 , NM_033138.3 , BC040354 , BC040354.1 , BC005006 , BC014035 , BC015839 , BC061926 , BF675655 , BF965359 , BQ010124 , BU683034 , BX538339 , BX648808 , BX649088 , NM_033138.4 |
| UniProt ID: | Q05682 |
| Synonyms: | CDM; H-CAD; HCAD; L-CAD; LCAD; NAG22 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)