

Product datasheet for **HP234238**

PAXBP1 Human qPCR Primer Pair (NM_016631)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	94104
Forward Sequence:	GATGGCACCAAATCTTGACTCCT
Reverse Sequence:	GTTCTCTGGCTTGCTACGACG
ACCN:	<u>NM_016631</u> , <u>NM_016631.1</u> , <u>NM_016631.2</u> , <u>NM_016631.3</u> , <u>BC005213</u> , <u>BC030539</u> , <u>BC062992</u> , <u>BC084542</u> , <u>BC156214</u> , <u>NM_016631.4</u>
UniProt ID:	<u>Q9Y5B6</u>
Synonyms:	BM020; C21orf66; FSAP105; GCFC; GCFC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)