

Product datasheet for **HP234252**

NARF Human qPCR Primer Pair (NM_012336)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26502
Forward Sequence:	GCTGTTCAACGAGGATGTGGAG
Reverse Sequence:	CAGCAAAGCGTAACACCACCTC
ACCN:	<u>NM_012336</u> , <u>NM_012336.1</u> , <u>NM_012336.3</u> , <u>BC016440</u> , <u>BC000438</u> , <u>BQ228233</u> , <u>NM_012336.4</u>
UniProt ID:	<u>Q9UHQ1</u>
Synonyms:	IOP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)