

Product datasheet for **HP234285**

SLC14A1 Human qPCR Primer Pair (NM_001128588)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6563
Forward Sequence:	TTGCTCAGCCAGGACAGGTCAT
Reverse Sequence:	CAGGGAGTAACAGCCACCAGAA
ACCN:	<u>NM_001128588</u> , <u>NM_001128588.1</u> , <u>NM_001128588.2</u> , <u>NM_001128588.3</u> , <u>BC050539</u> , <u>BC050539.1</u> , <u>BC040128</u> , <u>NM_001128588.4</u>
UniProt ID:	<u>Q13336</u>
Synonyms:	HsT1341; HUT11; JK; Jk(b); RACH1; RACH2; UT-B1; UT1; UTE
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)