

## Product datasheet for **HP234286**

### **NDRG2 Human qPCR Primer Pair (NM\_201535)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	57447
<b>Forward Sequence:</b>	CACTCTGTGGAGACACCATACG
<b>Reverse Sequence:</b>	CTGAAACAGTGGCTGGAAGCAAG
<b>ACCN:</b>	<a href="#">NM_201535</a> , <a href="#">NM_201535.1</a> , <a href="#">BC010458</a> , <a href="#">BC010458.1</a> , <a href="#">BC093038</a> , <a href="#">BC011240</a> , <a href="#">BC013209</a> , <a href="#">BF341160</a> , <a href="#">BG818977</a> , <a href="#">BI759006</a> , <a href="#">BM810623</a> , <a href="#">BX247987</a> , <a href="#">BX248031</a> , <a href="#">BX647091</a> , <a href="#">BX647748</a> , <a href="#">NM_201535.2</a>
<b>UniProt ID:</b>	<a href="#">Q9UN36</a>
<b>Synonyms:</b>	SYLD
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)