

Product datasheet for **HP234296**

PLD3 Human qPCR Primer Pair (NM_012268)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23646
Forward Sequence:	CGTCCTGCATACCAAGTTCTGG
Reverse Sequence:	GCTGCAGTTGTACATGACCACG
ACCN:	NM_012268 , NM_012268.1 , NM_012268.2 , NM_012268.3 , BC036327 , BC036327.1 , BC096820 , BC000553 , BP304863 , BQ640152 , NM_012268.4
UniProt ID:	Q8IV08
Synonyms:	AD19; HU-K4; HUK4; SCA46
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)