

## Product datasheet for **HP234321**

### SP140 Human qPCR Primer Pair (NM\_007237)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11262
Forward Sequence:	GCATGGCTGGAGCAGAATGAGA
Reverse Sequence:	CGTTCGCCTTCTTTTCACTGGAG
ACCN:	<a href="#">NM_007237</a> , <a href="#">NM_007237.1</a> , <a href="#">NM_007237.2</a> , <a href="#">NM_007237.3</a> , <a href="#">NM_007237.4</a> , <a href="#">BC038091</a> , <a href="#">BC054890</a> , <a href="#">BC070160</a> , <a href="#">BC105743</a> , <a href="#">BC105744</a> , <a href="#">BC105960</a> , <a href="#">NM_007237.5</a>
UniProt ID:	<a href="#">Q13342</a>
Synonyms:	LYSP100; LYSP100-A; LYSP100-B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)