

Product datasheet for **HP234325**

Nav1.5 (SCN5A) Human qPCR Primer Pair (NM_198056)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6331
Forward Sequence:	CAAGACCTGCTACCACATCGTG
Reverse Sequence:	GTCGGCATACTCAAGCAGAACC
ACCN:	<u>NM_198056</u> , <u>NM_198056.1</u> , <u>NM_198056.2</u> , <u>BC140813</u> , <u>BC051374</u> , <u>BC144621</u> , <u>BU845010</u> , <u>NM_198056.3</u>
UniProt ID:	<u>Q14524</u>
Synonyms:	CDCD2; CMD1E; CMPD2; HB1; HB2; HBBB; HH1; ICCD; IVF; LQT3; Nav1.5; PFHB1; SSS1; VF1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)