

Product datasheet for **HP234354**

PLA2R (PLA2R1) Human qPCR Primer Pair (NM_007366)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22925
Forward Sequence:	CCCAAGATTCCGTCTGGTACC
Reverse Sequence:	GATCACTGTGCAGCCAGCTACA
ACCN:	<u>NM_007366</u> , <u>NM_007366.1</u> , <u>NM_007366.2</u> , <u>NM_007366.3</u> , <u>NM_007366.4</u> , <u>BC140823</u> , <u>BC144631</u> , <u>BF725498</u> , <u>BF805380</u> , <u>BP351340</u> , <u>BQ025507</u> , <u>BQ637466</u> , <u>BX485732</u> , <u>NM_007366.5</u>
UniProt ID:	<u>Q13018</u>
Synonyms:	CLEC13C; PLA2-R; PLA2G1R; PLA2IR; PLA2R
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)