

## Product datasheet for **HP234364**

### Artemin (ARTN) Human qPCR Primer Pair (NM\_057160)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9048
Forward Sequence:	TGCTGAGCAGCGTCGCAGAGG
Reverse Sequence:	GCTCTTCCACTGCACCAGCGG
ACCN:	<u><a href="#">NM_057160</a></u> , <u><a href="#">NM_057160.1</a></u> , <u><a href="#">NM_057160.2</a></u> , <u><a href="#">BC062375</a></u> , <u><a href="#">BC062375.1</a></u>
UniProt ID:	<u><a href="#">Q5T4W7</a></u>
Synonyms:	artemin; ENOVIN; EVN; NBN; neublastin; neurotrophic factor; neurotrophic factor artemin; OTTHUMP00000009173; OTTHUMP00000009174
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)