

## Product datasheet for **HP234397**

### **EFTUD1 (EFL1) Human qPCR Primer Pair (NM\_024580)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	79631
<b>Forward Sequence:</b>	CTTCTGATGGGAAGGGAAGTGG
<b>Reverse Sequence:</b>	ATGGCAGGCTACACAGTGTTC
<b>ACCN:</b>	<u><a href="#">NM_024580</a></u> , <u><a href="#">NM_024580.1</a></u> , <u><a href="#">NM_024580.2</a></u> , <u><a href="#">NM_024580.3</a></u> , <u><a href="#">NM_024580.4</a></u> , <u><a href="#">NM_024580.5</a></u> , <u><a href="#">BC014628</a></u> , <u><a href="#">BC041565</a></u> , <u><a href="#">BC152910</a></u> , <u><a href="#">BC156743</a></u> , <u><a href="#">BG204990</a></u> , <u><a href="#">BX538332</a></u> , <u><a href="#">NM_024580.6</a></u>
<b>UniProt ID:</b>	<u><a href="#">Q7Z2Z2</a></u>
<b>Synonyms:</b>	EFTUD1; FAM42A; HsT19294; RIA1; SDS2
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)