

## Product datasheet for **HP234422**

### DIAPH3 Human qPCR Primer Pair (NM\_030932)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	81624
Forward Sequence:	GATGAAACACGGTTGGCAGAGTC
Reverse Sequence:	ACTGCTCAGGTTACATAAGTTGC
ACCN:	<a href="#">NM_030932</a> , <a href="#">NM_030932.1</a> , <a href="#">NM_030932.2</a> , <a href="#">NM_030932.3</a> , <a href="#">BC034952</a> , <a href="#">BC041395</a> , <a href="#">BC048963</a> , <a href="#">BC068504</a> , <a href="#">BF979672</a> , <a href="#">BG203073</a> , <a href="#">BX649186</a>
UniProt ID:	<a href="#">Q9NSV4</a>
Synonyms:	AN; AUNA1; DIA2; diap3; DRF3; mDia2; NSDAN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)