

Product datasheet for **HP234444**

ENTPD6 Human qPCR Primer Pair (NM_001247)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	955
Forward Sequence:	TGCTACTCCAGGAGTTCGGCTT
Reverse Sequence:	GCCACTATGAGGCTGGACTCTT
ACCN:	NM_001247 , NM_001247.1 , NM_001247.2 , NM_001247.3 , NM_001247.4 , BC025980 , BC025980.2 , BM685060 , BM982256 , NM_001247.5
UniProt ID:	O75354
Synonyms:	CD39L2; dj738P15.3; IL-6SAG; IL6ST2; NTPDase-6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)