

Product datasheet for **HP234447**

ZNF451 Human qPCR Primer Pair (NM_015555)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26036
Forward Sequence:	CTCTCTGTGCAAAGATGTTCCCT
Reverse Sequence:	CAGCTACAGGTCCAGCATTTCCG
ACCN:	<u>NM_015555</u> , <u>NM_015555.1</u> , <u>NM_015555.2</u> , <u>BC021712</u> , <u>BC042450</u> , <u>BC058853</u> , <u>BC098382</u> , <u>BM559241</u> , <u>BX648081</u>
UniProt ID:	<u>Q9Y4E5</u>
Synonyms:	COASTER; dj41711.1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)