

Product datasheet for **HP234510**

KCNK2 Human qPCR Primer Pair (NM_014217)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3776
Forward Sequence:	CTGCTGTCCTGAGCATGATTGG
Reverse Sequence:	TGTGACGTTGGCTGTCCACTCA
ACCN:	<u>NM_014217</u> , <u>NM_014217.1</u> , <u>NM_014217.2</u> , <u>NM_014217.3</u> , <u>BC069462</u> , <u>BC069462.1</u> , <u>BC101693</u> , <u>BC101695</u> , <u>BC143586</u> , <u>NM_014217.4</u>
UniProt ID:	<u>O95069</u>
Synonyms:	hTREK-1c; hTREK-1e; K2p2.1; TPKC1; TREK; TREK-1; TREK1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)