

Product datasheet for **HP234605**

FAM38A (PIEZO1) Human qPCR Primer Pair (NM_001142864)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9780
Forward Sequence:	CCTGGAGAAGACTGACGGCTAC
Reverse Sequence:	ATGCTCCTTGGATGGTGAGTCC
ACCN:	<u>BC008073</u> , <u>NM_001142864</u> , <u>NM_001142864.1</u> , <u>NM_001142864.2</u> , <u>NM_001142864.3</u> , <u>BC008073.1</u> , <u>BC016983</u> , <u>BC141877</u> , <u>BC150271</u> , <u>BE066337</u> , <u>BU633404</u> , <u>NM_001142864.4</u>
UniProt ID:	<u>Q92508</u>
Synonyms:	DHS; FAM38A; LMPH3; LMPHM6; Mib
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)