

Product datasheet for **HP234663**

CAP1 Human qPCR Primer Pair (NM_006367)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10487
Forward Sequence:	CAGTCTCTACCACTTCAGGCTC
Reverse Sequence:	ACTGGACCACTCTGAGCCTCA
ACCN:	NM_006367 , NM_006367.1 , NM_006367.2 , NM_006367.3 , BC013963 , BC013963.2 , BC095440 , BC020367 , BI460362 , BM772471 , BP391532 , BT007152 , BX537801
UniProt ID:	Q01518
Synonyms:	CAP; CAP1-PEN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)