

Product datasheet for **HP234664**

BIN1 Human qPCR Primer Pair (NM_139343)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	274
Forward Sequence:	CGTCAACACGTTCCAGAGCATC
Reverse Sequence:	CTTGACCGTGAAGGTGTTGCTC
ACCN:	<u>NM_139343</u> , <u>NM_139343.1</u> , <u>NM_139343.2</u> , <u>BC004101</u> , <u>BE047931</u> , <u>BM711160</u> , <u>BM742027</u> , <u>BM805510</u> , <u>BT006865</u> , <u>NM_139343.3</u>
UniProt ID:	<u>O00499</u>
Synonyms:	AMPH2; AMPHL; CNM2; SH3P9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)