

## Product datasheet for **HP234777**

### Factor VIII (F8) Human qPCR Primer Pair (NM\_019863)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2157
Forward Sequence:	GCCATTGGGAATGGAGAGTAAAG
Reverse Sequence:	CCTGAGTAGTTACTCCTGTGAC
ACCN:	<a href="#">NM_019863</a> , <a href="#">NM_019863.1</a> , <a href="#">NM_019863.2</a> , <a href="#">BC022513</a> , <a href="#">BC064380</a> , <a href="#">BC098389</a> , <a href="#">BC111967</a> , <a href="#">BC111969</a> , <a href="#">BC150305</a> , <a href="#">BC166700</a> , <a href="#">BI547023</a>
UniProt ID:	<a href="#">P00451</a>
Synonyms:	AHF; DXS1253E; F8B; F8C; FVIII; HEMA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)